

# Simon<sup>dr</sup>

D I G I T A L   R A D I O G R A P H Y



VETSTANDARD  
PREMIUM X-RAY GENERATION  
BY QUANTUM MEDICAL IMAGING

DRV  
HIGH-RESOLUTION DIGITAL X-RAY  
BY SIMONDR, INC.

From a leader in analog radiology  
for more than 50 years, SimonDR is  
lighting the way to the future of  
veterinary digital radiography.

# VETSTANDARD

PREMIUM X-RAY GENERATION  
BY QUANTUM MEDICAL IMAGING  
WITH DIGITAL X-RAY INTEGRATION BY SIMONDR, INC.

A COMPLETE NO-COMPROMISE  
**DIRECT RADIOGRAPHY SYSTEM**  
THAT'S EASY TO USE



## FEATURES AND OPTIONS

### **VetStandard Features Include:**

*SimonDR DRV high-resolution CCD based direct digital x-ray system.*  
*3.5 line pairs per millimeter.*  
*13 million pixel sensor.*  
*Made and supported in the U.S.A.*  
32kw - 25ma - 500ma high-frequency standard power.  
Higher ma values mean faster exposure times.  
200,000 heat-unit tube with a standard fractional focal spot (.6mm / 1.5mm).  
60 inch four-way floating table top with 30 inches of longitudinal float to accommodate larger animals.  
Included "Uro Track" for easy liquid containment.  
Patient Tie-downs.  
Recessed foot-control operated floating table-top with premium fail-safe brakes.  
650 pound patient capacity tabletop.  
Laser guided collimator for ultimate accuracy.  
Pre-set anatomical exposure factors capable.

### **VetStandard Warranty:**

2 year receptor and standard 5 year limited x-ray generation warranty.

### **Options Include:**

Stored energy gives you the ability to use high-frequency where incoming power is limited.  
Shorter float for very small areas.  
Tubestand-mounted preview monitor.

**Simon***dr*

DIGITAL RADIOGRAPHY

800-835-3852

3515 MARMENCO COURT  
BALTIMORE, MARYLAND 21230