




Austin's X-Ray
Service, Inc.

Chiropractic C-Stand



The Chiropractic C-Stand, from Austin's X-Ray Service, is designed specifically for the Chiropractic community. The C-arm design maintains constant alignment between the x-ray tube and image receptor in both vertical and horizontal positions. This flexibility allows x-rays to be taken on patients in standing, sitting, or laying positions.

The Chiropractic C-Stand is also designed to facilitate the use of any digital receptor or bucky. This feature makes the C-Stand a perfect solution for chiropractors who are converting to digital x-ray, or who are anticipating a move to digital x-ray in the future.

The Chiropractic C-Stand can be purchased as a package with any of our digital or analog solutions, or as an individual product for use with your current equipment.

The C-arm design comes standard with a 40 inch source-to-image distance (SID), with an option to upgrade to a model that can be used with either a 40 or 72 inch SID.

AUSTIN'S X-RAY SERVICE
5266 COUNTY ROAD 2
SWANTON, OH 43558

Phone: 419-825-3834
Email: info@austinsxrayservice.com
Website: www.austinsxrayservice.com